




AUTOCLAVE-VERTICAL

STERILIZER FOR DRY STERILIZATION (VERTICAL)

Outer Body made of thick Stainless Steel/Mild Steel sheet finished with oven baked enamel. Inner chamber will be thick Stainless Steel. The chamber is absolutely leakproof and it can be operated at any selective point in between 5 to 20 psi (pounds/sq. inch). The instrument is equipped with pressure gauge, pressure release valve, safety valve, water level indicator and perforated S.S. basket complete with cord and plug. The lid is made of thick Stainless Steel Sheet with pedal lifting arrangement and radial locking devices. To work on 220/240V, AC Single phase supply.



Inner Chamber Dia × Depth

- a) 250×450 MM
- b) 300×500 mm
- c) 350×550 mm
- d) 400×600 mm
- e) 450×600 mm
- f) 550×750 mm

Load

- 1.5 kw
- 2.0 kw
- 3.0 kw
- 4.0 kw
- 4.0 kw
- 5.0 kw

Optional Accessories :

- 1. Automatic low water level cut-off device
- 2. Pressure control switch
- 3. Automatic 0-120 minute timer (Imported)
- 4. Digital Temperature Indicator

Mktd by : **AIM SCIENTIFIC**

e-mail : info@aimscientific.com
Website : www.aimscientific.com