



## HOT PLATE RECTANGULAR

### SPECIFICATIONS:

- Used for continuous heating upto 350°C
- Top Rectangular Cast Iron plate fitted with Insulated heating elements.
- Temperature is controlled by Energy Regulator or thermostat.
- Outer Body made of Mild Steel painted with white stoving enamel.
- Complete with pilot Indicators, cord & plug. To operate on 220/230 V 1 Phase, AC Supply.

S.No.	Available Sizes (L×B×H)	Load
a)	30×25×15 cm.	(1.2 Kw.)
b)	40×25×15 cm.	(1.5 Kw.)
c)	45×30×15 cm.	(2.0 Kw.)
d)	60×45×15 cm.	(3.0 Kw.)

### Optional Accessories :

Digital Temperature Controller cum Indicator

### SPECIFICATIONS:

- Used for continuous heating upto 350°C
- Round Cast Iron Plate, 20.5 cm dia. fitted with Insulated heating elements. (1.25 Kw)
- Temperature is controlled by 3 heat control switch or energy regulator/thermostat.
- Outer Body made of Mild Steel painted with white stoving enamel.
- Complete with pilot Indicators, cord & plug. To operate on 220/230 V 1 Phase, AC Supply.

Special Sizes can be made as per user requirements  
Right reserved to change specifications.

## LABORATORY HOT PLATE ROUND

