



HEATING MANTLE

SPECIFICATIONS:

- Used for continuous heating upto 350°C of inflammable liquids.
- Electric heating net is jointless, hand knitted from glass yarn.
- The Body of heating mantle is spun one piece made from Mild Steel painted with white stoving enamel.
- Heating Mantles are logged with special mineral wool.
- Temperature is controlled by energy regulator/ thermostat.
- Complete with pilot Indicators, cord & plug.

S.No.	Available Sizes (L×B×H)	Load
a)	50 & 100 ml.	(80 watt.)
b)	250 & 500 ml.	(200 watt.)
c)	01 Ltr.	(300 watt.)
d)	02 Ltr.	(450 watt.)
e)	05 Ltr.	(600 watt.)
f)	10 Ltr.	(2×600 watt.)
g)	20 Ltr.	(2×650 watt.)



MAGNETIC STIRRER without/with hot plate

SPECIFICATIONS:

- Compact Stirring device utilizing a rotating field of magnetic induce variable speed stirring action with either closed or open vessel.
- Stirring is accomplished by means of small rotor sealed in teflon which is placed in the liquid capable of rotating by magnetic field applied from below.
- Speed is Variable upto 1200 RPM and controlled by Speed regulator.
- Stainless Steel Top.
- Complete with pilot Indicator, cord & plug.

S.No.	Available Sizes	Flask Capacity
a)	Capacity	(1.0 Ltr.)
b)	Capacity	(2.0 Ltr.)
c)	Capacity	(5.0 Ltr.)
d)	Capacity	(10.0 Ltr.)

Optional Model : Magnetic Stirrer with Hot Plate fitted with energy regulator to control Temperature.